

CLOVER⁺ A1cTM Self

HbA1c Monitoring System

Fully
Automatic

4 μ l

5 min

- Fully Automatic
- Boronate Affinity
- Spectrophotometry
- 4 μ l Blood Sample
- Test Result in 5 Minutes



INFOPIA Co., Ltd.

132, Anyangcheondong-ro, Dongan-gu, Anyang-si, Gyeonggi-do, Republic of Korea
www.infopia21.com  www.facebook.com/infopia21




IMM21-020 Rev.2014-01-04

CLOVER[®] A1c[™] Self

Accuracy

1. Comparison method

CLOVER A1c[™] Self analyzer vs Tosoh G8, N=120




| Sample type | Regression line | slope | y-interceptor | R2 |
|-----------------|-----------------|--------|---------------|--------|
| capillary blood | y=0.9835x-0.004 | 0.9835 | -0.004 | 0.9803 |

Table1. Summary of regression line

2. System bias plot (Total)

Bias plot (N=120, Limit: within bias $\pm 7\%$)



| Within $\pm 5\%$ | Within $\pm 7\%$ | Within $\pm 10\%$ |
|------------------|------------------|-------------------|
| 99/120(82.5%) | 114/120(95.0%) | 119/120(99.2%) |

Table2. System accuracy results of capillary for HbA1c levels (mean bias, %)

Specification

| | |
|-------------------------|--|
| Sample type | Capillary whole blood, Venous blood with anticoagulant |
| Sample volume | 4 μ l |
| Test range | 4.0~14.0% |
| Reading time | 5 minutes |
| Memory capacity | 200 test results |
| Wavelength | 415 nm (normal) |
| Power required | DC 9V-1.5 A |
| Dimensions / Weight | 200 x 198 x 139(mm) / 1.4kg |
| Storage temperature | 10~35 $^{\circ}$ C (50~95 $^{\circ}$ F) |
| Operating temperature | 17~32 $^{\circ}$ C (63~90 $^{\circ}$ F) |
| Relative humidity range | 10%~90% |

Precision

1. Acceptance criteria of CV

Total imprecision (CV) must not be statistically significantly 4%

2. Test result

| | Normal | Abnormal |
|-----------|--------|----------|
| Mean, (%) | 5.3 | 10.0 |
| SD | 0.14 | 0.19 |
| CV, (%) | 2.7 | 1.9 |

Table3. Mean values, Standard deviation (SD), coefficient variation (CV) of each sample

3. Conclusion

All CV at each HbA1c level are in the acceptance criteria.

Test Procedure




01 Open the CLOVER A1c[™] Self Analyzer.



02 Insert the Cartridge.



03 Shake the Reagent Pack 5~6 times




04 Place the blood sample (4 μ l) at the sample Collecting Area.




05 Insert the Reagent Pack into the Cartridge in the Analyzer.



06 Close the cover of Analyzer.



07 The test starts automatically by closing the cover, and the test result time countdown is displayed.



08 After 5 minutes, when the test is over, the HbA1c% is displayed.

Full Package

CloverA1c[™] Self analyzer, Instruction for Use, Exclusive power adapter, Fan Filters, Daily check cartridge, Monthly check cartridge

**Options : Blood Glucose Module, PC connection cable, Barcode Scanner, Thermal Printer